



Stellant® Sx CT Injection System

Your Intelligent, Total
Fluid Delivery System

Intelligent, Intuitive, Innovative

For four decades, MEDRAD has been changing the way medicine is practised by developing some of the most technologically-advanced medical devices in health care. Now we continue our tradition of contrast delivery innovation with the Stellant Sx Injection System. Stellant Sx is the result of listening to and meeting the needs of health care professionals in the most demanding CT environments.

MEDRAD's Stellant Sx is the Right Choice for the CT Environment of Today and Tomorrow



Prime Button

- Pushes a predetermined amount of contrast into the tubing to displace air.

Remote Check For Air

- Confirm that no air is in fluid path at the injector head to reduce risk of an air injection.

Remote Arming

- Arm the injection system at the injector head.

Improved Syringe Design Saves Time and Increases Throughput Capability

Quick and Easy Syringe Installation

- Simply snap the syringe on in any orientation.
- Plunger advances automatically if contrast syringe is empty.
- Just twist to remove the syringe following the procedure – plunger will retract automatically into position for the next case.


Integral Auto-load

- Simply connect the fluid path, and the syringe is automatically loaded with the volume of contrast required for the protocol.



Priming Tube

- Mini-reservoir captures contrast overflow while priming patient tubing.



“The Stellant Sx Single Xpanded technology has really helped improve the quality of CT Angiography studies that we do.”

The New Stellant User Interface Provides Many Powerful Features

The Stellant Interface was designed from the ground up to be the most intuitive and user-friendly system available.



Colour Touch Screen

- State-of-the-art colour touch screen makes all system operations more intuitive and easier to view.
- Easy-to-read, anti-glare screen.
- Move easily from page to page without scrolling.
- Stores and recalls up to 32 different protocols.



Pressure Monitor Graph


- Displays the pressure being developed in the syringe during an injection.
- Allows operator to reduce patient risk by providing real-time visual graph of pressure.

Flow Profile

- Provides a real-time graphic display of the entire injection profile. You will always know at a glance what has been programmed, and where the injection system is in the delivery of the protocol.

Stop-Watch Feature

- Timer starts when the injection is initiated. This feature is especially useful to time the start of a delayed phase scan.



Scanning Delay

- Enables pre-set or injection delay to maximise synchronisation with the scanner.

Multi-phase Option

- Up to six phases and Hold & Pause options permit full capability for today's complicated scan protocols.

Protocol Lock

- Current protocol can be locked at the display to enable remote arming and injections from the patient's side to enhance safety.

Performance. For life.

Since 1964, MEDRAD has been developing products that clinicians worldwide rely on to provide outstanding performance every time they're used. Performance that helps physicians diagnose patients more accurately and, ultimately, deliver quality patient care. You will find that same level of performance – along with the high level of customer support you've come to expect from MEDRAD – with the new Stellant Sx Computed Tomography Injection System.



System Specifications

Injection Specifications

Flow Rate (range & increments)	0.1 to 10 mL/sec in 0.1 mL increments
Volume (range & increments)	1 mL to syringe capacity in 1 mL increments
Programmable Pressure Limit psi/kPa	200 mL syringe: 325 psi, 2241 kPa
Scanning Delay	0-300 seconds (5 minutes) in 1 second increments
Pause	1-900 seconds (15 minutes) in 1 second increments
Hold	Maximum HOLD time is 20 minutes
Syringes (Volume capacity)	200 mL sterile disposable syringe
Maximum Number of Phases	6
Maximum Number of Protocols	32

Miscellaneous Specifications

Electrical Requirements (VAC/Hz)	100-240 VAC, 50/60 Hz
Syringe Heater Range	95 degrees F +/-9 degrees 35 degrees C +/-5 degrees
Head Mounting Options	Overhead Counterpoise System or Floor Pedestal (Floor Pedestal includes an Integral IV Pole)

Mechanical Specifications

	Height	Width	Depth	Weight
Single Injector Head	6.1" (15,5 cm)	8.6" (21,8 cm)	14.5" (36,8 cm)	12 lbs (5,4 kg) <i>(without syringe)</i>
Display Control Unit with stand	13.5" (34,3 cm)	12.5" (31,9 cm)	8.9" (22,5 cm) <i>Screen only = 3.4" (8,6 cm)</i>	8.2 lbs (3,7 kg)
Base Unit	11.5" (29,2 cm)	11.0" (27,9 cm)	8.8" (22,2 cm)	13.6 lbs (6,2 kg)

Ordering Information

Systems listed include all standard cabling

Catalogue Number: SCT 121

Catalogue Number: SCT 122

Stellant Sx Injection System
with Pedestal Head Mount

Stellant Sx Injection System
with Counterpoise System

Backed by MEDRAD reliability, support and service

- Legendary product reliability
- Unmatched after-sales support
- Knowledgeable and reliable technical service
- Ask about MEDRAD's Predictive Maintenance Programs and Extended Warranties

MEDRAD reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized MEDRAD representative for the most current information.

MEDRAD and Stellant are federally registered trademarks and MEDRAD Radiology and Performance. For life. are trademarks of MEDRAD, INC. U.S.A.

MEDRAD, INC.
One Medrad Drive
Indianola, PA 15051-0780 USA
Telephone: +1-412-767-2400
Fax: +1-412-767-4120

Customer Service/Orders
+1-800-MEDRAD-1
+1-800-633-7231

Customer Service Fax
+1-412-767-4120

www.medrad.com
Email address for more information:
info@medrad.com