



## Flexibility, Ease-of-Care and Increased Throughput with Certo® MR Wireless Network



**Certo MR Wireless Network provides reliable communication throughout the MR suite, enabling flexibility, ease-of-care, and increased throughput.**

### **Flexibility, Ease-of-Care, and Increased Throughput**


- Certo MR Wireless Network provides greater flexibility by allowing MEDRAD wireless-enabled equipment to move freely both in the MRI scan room and in the control room.
- Clinicians can focus on patient care while efficiently managing patients through the Veris® MR Patient Monitoring System.
- Certo allows Veris to freely travel with the patient from the prep room to the scanner room and then to the recovery area while maintaining communication with the remote display on the network.

### **Performance and Safety**

- Fiber optic connection between scanner room and control room doesn't interfere with imaging and provides a secure path for transmitting your patient's vital sign information.

*MEDRAD has a philosophy of taking care of the patient,  
to better the outcome for the patient.\*\**

**Certo® MR Wireless Network**



**Innovative Solutions, Reliable Technology**

Certo offers an innovative solution to wireless communication in the MR environment while utilizing the well-established 802.11 b/g wireless networking standard featuring a 2.4 GHz Radio Frequency (RF) interface.

Contact your MEDRAD Sales Representative for more information at 1-800-MEDRAD-1

**Backed by MEDRAD reliability, support and service**

- Legendary product reliability
- Unmatched after-sales support
- Knowledgeable and reliable technical service
- Ask about MEDRAD's Predictive Maintenance Programs and Extended Warranties

MEDRAD reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized MEDRAD representative for the most current information.

MEDRAD, Certo and Veris are federally registered trademarks, and MEDRAD Radiology and Performance. For life., are trademarks of MEDRAD, INC. U.S.A.

\*\*Based upon 2003 - 2004 FSC (Flaherty, Sabol & Carroll) Communications research.

**MEDRAD, INC.**  
 One Medrad Drive  
 Indianola, PA 15051-0780 USA  
 Telephone: +1-412-767-2400  
 Fax: +1-412-767-4120

**www.medrad.com**  
 Email address for more information:  
 info@medrad.com

**Customer Service/Orders**  
 +1-800-MEDRAD-1  
 +1-800-633-7231

**Customer Service Fax**  
 +1-412-767-4120

