

MEDRAD CT

Stellant® Interface Module
Simplify your scan.

*"That's what I know about MEDRAD,
it works."*



Synchronized communication between your Stellant CT Injection System and the CT scanner from one location is within your grasp.



- Request a CT scanner start or stop remotely with Stellant or request an injection system start or stop with the CT scanner
- Automatically request that the CT scanner begin its scanning sequence for proper timing with a Stellant programmed scan delay
- Send a request to simultaneously stop both your Stellant and the CT scanner in emergency situations

MEDRAD'S Stellant Imaging System Interface (ISI) Module makes it simple for your Stellant and the scanner to communicate.

Stellant ISI delivers precisely-timed synchronized communication. It enables the start of your Stellant Injector and a request to start the CT scanner from a single location - no more back and forth from room to room.

Simplify your scan with Stellant ISI.

Stellant ISI is currently supported by Philips Brilliance and Mx8000 Scanners containing SAS and requires appropriate cable from scanner manufacturer.

MEDRAD®
Performance. For life.™

MEDRAD reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized MEDRAD representative for the most current information.

MEDRAD, and Stellant are federally registered trademarks, and Performance. For life. is a trademark of MEDRAD, INC. U.S.A.

Philips is a federally registered trademark, and Brilliance is a trademark of Philips Medical Systems, U.S.A. Used by permission.



© 2008 MEDRAD, INC. All Rights Reserved.

STL ISI NA 202690 Rev B



The Malcolm Baldrige National Quality Award is our nation's premier award for performance excellence and quality achievement. Established in 1987 and named after the former Secretary of Commerce, the Baldrige Award recognizes exemplary achievements in seven areas: leadership, strategic planning, customer and market focus, measurement, analysis and knowledge management, human resource focus, process management and business results. The U.S. Commerce Department's National Institute of Standards and Technology (NIST) manages the Baldrige National Quality Program in close cooperation with the private sector. The award is traditionally presented by the President of the United States in a special Washington, D.C. ceremony.

MEDRAD, INC.
One Medrad Drive
Indianola, PA 15051-0780 USA
412-767-2400

Customer Service/Orders
1-800-MEDRAD-1
(1-800-633-7231)

Customer Service FAX
(412) 767-4120

International FAX
(412) 767-4128

For Service Repairs
1-800-MEDRAD-7
(1-800-633-7237)

MEDRAD Service FAX
(412) 767-4126

www.medrad.com/stellant
Email address for more
information:
info@medrad.com

MEDRAD®
Performance. For life.™