

Intego™ PET Infusion System Specifications

Vial Compatibility (March 2010)

Vial Shield and Vial Specifications							
Vial Shields			Vials				
MEDRAD P/N	Collar Color	Weight	Volume Range	ISO 8362-1:2003 ISO 8362-4:2003 Size Specification	Body Diameter (mm)	Glass Height (mm)	Cap Diameter (mm)
3017469	Silver	18 lbs (8.6 kg)	30-35 ml	30H	36 ± 0.5	62.8 ± 0.7	20
3025071	Blue	18 lbs (8.6 kg)	20-25 ml	20R	30 ± 0.25	55 ± 0.7	20
3027139	Gold	18 lbs (8.6 kg)	20-25 ml	20H	32 ± 0.45	58 ± 0.6	20
3025939	Green	18 lbs (8.6 kg)	10-15 ml	10R	24 ± 0.2	45 ± 0.5	20
3025940	Brown	18 lbs (8.6 kg)	10-15 ml	15R	24 ± 0.25	60 ± 0.5	20
3027305	Black	18 lbs (8.6 kg)	10-15 ml	10H	25 ± 0.25 25.4 ± 0.4 25.5 ± 0.25	53.5 ± 0.5	20

