

*Medrad announces the*  
**NEW Adjustable  
T-Rail Table  
Bracket**

The new adjustable T-Rail Table Bracket accommodates most standard x-ray tables. The table bracket is compatible with Medrad *Mark V ProVis™*, *Mark V Plus™*, *Mark V™* and *Mark IV™* injector heads. The new table bracket adjusts easily with the turn of a knob to a convenient height for the user. A simple press of a lever allows for quick, easy installation onto any standard x-ray table.

- *Height adjustable*
- *Release with push of a lever*
- *Easy four-step installation*
- *One hand adjustment*



# A N N O U N C E M E N T

## PRODUCT INFORMATION

CATALOG #	DESCRIPTION	LIST PRICE
KMA 350	Universal Adjustable T-Rail Table Bracket with One Year Warranty	\$ 1,495.00

FEATURE	BENEFIT
Height adjustable: 11.75" (30cm) - 17" (43cm).	Adapts to user, procedure and room requirements. Improved safety (reduces potential of collision with imaging equipment).
Release with push of a lever.	Convenient to remove T-Rail bracket from Table Rail. Ability to change quickly.
Easy four-step installation.	No time wasted in set-up. Helps maximize efficiency.
Accommodates rails from .25" (6,4mm) to .5" (12,7mm) thick and .875" (22,2mm) to 1.125" (28,8mm) high.	Flexibility—can be used on all major imaging system tables.
Accepts Medrad <i>Mark IV</i> , <i>Mark V</i> , <i>Mark V Plus</i> , and <i>Mark V ProVis</i> injector heads.	Accommodates upgrades.
Stainless steel construction.	Strong and durable for reliable performance. Easy to clean.
No tools required for height adjustment.	No time wasted in setting up for a new user or procedure.

Medrad reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Medrad representative for the most current information.

## **MEDRAD**<sup>®</sup>

**Medrad, Inc.**  
**Corporate Office**

One Medrad Drive  
Indianola, PA 15051-0780 U.S.A.  
Telephone: (412) 767-2400  
FAX: (412) 767-4128

Medrad is a federally registered trademark and Mark V ProVis, Mark V Plus, Mark V and Mark IV are trademarks of Medrad, Inc. U.S.A.