

Configuration File Information MEDRAD 1.5T Breast Array

9x Software

The Configuration File for the Medrad 1.5 T Breast Array Coil is located in the system's Configuration File Manager. Please locate the following file in the Available Coil List:

BREASTPA This file turns on both sides of the coil

Click and Add this file name to the site coil list. A system re-boot is necessary following the addition of this file. The coil will then be available in the phased array coil selections.

Please call 1800 MRCOILS with any questions

Please note:

The ability to scan with the left only or right only side of the Medrad Breast Array Coil is NOT available until the Excite Platform.